



SPECIFICATION SHEET

Product Name : Sodium hydroxide 0.25 mol/l (0.25 N) , Volumetric solution
Product Number : 969.21V
Brand : ISOLAB
CAS Number : -
Formula : -
Formula Weight : -

| TEST NAME | SPECIFICATION | UNIT |
|-----------------------------------|-------------------------|-------|
| Amount of substance concentration | 0.24875 - 0.25125 mol/l | mol/l |
| Traceability | NIST | |
| Indicator | phenolphthalein | |

Quality Management
Isolab Laborgerate GmbH

