



SPECIFICATION SHEET

Product Name : tert-Butyl methyl ether, Extra pure
Product Number : 908.022
Brand : ISOLAB
CAS Number : 1634-04-4
Formula : C₅H₁₂O
Formula Weight : 88.15 g/mol

TEST NAME	SPECIFICATION	UNIT
Purity (GC)	≥ 99.0 %	%
Identity (IR)	passes test	
Free acid (as CH ₃ COOH)	≤ 0.002 %	%
Density (d 20°C/ 4°C)	0.740 - 0.742	gr/ml
Peroxide (as H ₂ O ₂)	≤ 0.005 %	%
Evaporation residue	≤ 0.005 %	%
Water	≤ 0.05 %	%

Quality Management
Isolab Laborgerate GmbH

