



SPECIFICATION SHEET

Product Name : Mercury(II) sulfate, Extra pure
Product Number : 947.032
Brand : ISOLAB
CAS Number : 7783-35-9
Formula : HgSO_4
Formula Weight : 296.65 g/mol

TEST NAME	SPECIFICATION	UNIT
Assay (complexometric)	$\geq 99.0 \%$	%
Identity	passes test	
Chloride (Cl)	$\leq 0.003 \%$	%
Fe (Iron)	$\leq 0.005 \%$	%
Mercury (I) (as Hg)	$\leq 0.15 \%$	%

Quality Management
Isolab Laborgerate GmbH

